

## BOB STANDING-FRAME



BOB es un standing-frame que permite ir desde la posición horizontal hasta la posición vertical a través de mando eléctrico, permitiendo asumir la posición supino. Con sección de tronco y pélvis ajustable en altura a través de columna central, este standing-frame permite vários ajustes independientes al nível de los soportes de rodillas y pies. Soporte de cabeza 100% ajustable. Cinturones de doble contención para mayor seguridad. Tablero (opcional) ajustable en altura y inclinación. Posibilidad de colocación de diversos accesorios.

*BOB is a standing frame that allows to go from the horizontal to the vertical position through electric control, allowing to assume the supine position. With trunk and pelvis section ajustable in height through central column, this standing frame also allows for several independent adjustments in knee and foot supports. 100% adjustable headrest. Double containment belts for greater security. Board (optional) ajustable in height and inclination. Several accessories available.*

**ORTHOSXXI, Unipessoal Lda.**

Rua Santa Leocádia nº2735 // 4809-012 Sta. Leocádia de Briteiros // Guimarães // Portugal

Phone + 351 253 470 276 // Mobile +351 936 199 005 // Email [margarida.silva@orthosxxi.com](mailto:margarida.silva@orthosxxi.com) // Website [www.orthosxxi.com](http://www.orthosxxi.com)


**Características / Characteristics**



**Permite varias posiciones**  
*Allows several positions*



**Cinturones de doble contención**  
*Double containment belts*



**Totalmente ajustable**  
*Fully adjustable*



**Actuador electrico**  
*Electrical actuator*

BOB standing-frame			
Marca/ Brand: <b>Orthos XXI</b>			
Longitud base <i>Base length</i>	1100mm	Soporte lateral tronco <i>lateral trunk support</i>	170x110mm
Ancho base <i>Base width</i>	600mm	Soporte lateral pélvis <i>Pelvic lateral support</i>	170x110mm
Soporte pélvico <i>Pelvic support</i>	380x180mm	Tablero (opcional) <i>Board (optional)</i>	650x550mm
Soporte de tronco <i>Trunk support</i>	380x180mm	Peso máximo usuario <i>Maximum user weight</i>	75Kg

Garantía/ Warranty
2 años contra defectos de fabricación <i>2 years against manufacturing defects</i>
Tiempo estimado de entrega / Delivery time
Consultar/ Consult

**ORTHOSXXI, Unipessoal Lda.**

Rua Santa Leocádia nº2735 // 4809-012 Sta. Leocádia de Briteiros // Guimarães // Portugal

Phone + 351 253 470 276 // Mobile +351 936 199 005 // Email [margarida.silva@orthosxxi.com](mailto:margarida.silva@orthosxxi.com) // Website [www.orthosxxi.com](http://www.orthosxxi.com)