

AZTECA COMPACT



AZTECA COMPACT es una silla eléctrica compacta y fácil de transportar. Concebida para uso en espacios estrechos, en interiores y pisos más regulares. Con estructura plegable y respaldo partido, para fácil almacenaje. Apoya brazos abatibles. Reposapiés removibles y ajustables en altura. Tapizados negros. Motores de 100w y baterías de 26Ah. Ruedas traseras y delanteras macizas de 200mm. Velocidad máxima hacia delante: 6Km/h. Peso máximo del usuario: 125Kg.

AZTECA COMPACT, is a electric wheelchair with compact dimensions and easy to transport. Conceived for use in narrow spaces, interiores and regular floors. With folding structure and folding back, for easy storage. Folding armrests. Removable footrests, height adjustable. Black upholstery. 100w motors and 26Ah batteries. 200mm rear and front solid wheels. Maximum speed forward: 6Km/h. Maximum user weight: 125Kg.

ORTHOSXXI, Unipessoal Lda.

Rua Santa Leocádia nº2735 // 4809-012 Sta. Leocádia de Briteiros // Guimarães // Portugal

Phone + 351 253 470 276 // Mobile +351 936 199 005 // Email margarida.silva@orthosxxi.com // Website www.orthosxxi.com

Características / Characteristics



Estructura compacta
Compact structure



Respaldo partido
Folding back



Reposapiés desmontables
Removable footrests



Joystick compacto
Compact joystick

Azteca Compact			
Marca/ Brand: Orthos XXI			
Ruedas traseras <i>Rear wheels</i>	200mm	Ancho asiento <i>Seat width</i>	40,43,46cm
Ruedas delanteras <i>Front wheels</i>	200mm	Profundidad asiento <i>Seat depth</i>	42cm
Estructura <i>Structure</i>	Negro <i>Black</i>	Motores <i>Motors</i>	100W
Tapizados <i>Upholstery</i>	Negro <i>Black</i>	Baterías	12V 26Ah
Velocidad máxima hacia delante <i>Maximum speed forward</i>	6Km/h	Autonomía <i>Autonomy</i>	20km
Velocidad máxima hacia atrás <i>Maximum reverse speed</i>	3Km/h	Peso max.usuario <i>Maximum user weight</i>	125Kg

Garantía/ Warranty
2 años contra defectos de fabricación <i>2 years against manufacturing defects</i>
Tiempo estimado de entrega / Delivery time
Consultar/ Consult

ORTHOSXXI, Unipessoal Lda.

Rua Santa Leocádia n°2735 // 4809-012 Sta. Leocádia de Briteiros // Guimarães // Portugal

Phone + 351 253 470 276 // Mobile +351 936 199 005 // Email margarida.silva@orthosxxi.com // Website www.orthosxxi.com